

LEUKOMED® CONTROL

DEAR PATIENT,

your wound is securely protected by **Leukomed® Control**.

This dressing is designed by our Leukoplast® experts to provide ideal healing conditions for your skin and wound. Here are a few tips to ensure that **Leukomed® Control** keeps protecting your skin and wound.



Leukomed® Control is waterproof to protect your skin. It allows you to take quick showers. We recommend that you avoid prolonged soaking or bathing, though.



Leukomed® Control feels like a second skin, smoothly conforming to your body's movements. Nevertheless it is recommended to refrain from physical activities that cause intense sweating or extreme muscle flexing until the skin has fully regained its strength.



Leukomed® Control is highly transparent, so that you and your doctor can easily inspect the progress of your wound healing. We recommend that you do not put any creams, lotions or fluid medication on or underneath the dressing.



Your doctor will tell you when **Leukomed® Control** can be removed. If you experience exceeding pain, swelling, redness or warming of the wound surrounding or if the wound opens up, please contact your doctor immediately.



Leukomed Control

Sterilen, polpropusten vodoodporen poliuretanski film s hidropolimerno blazinico za vpijanje izločka kirurških in drugih akutnih ran.

Omogoča hiter pregled nad stanjem rane in izločkom, ter vzpodbuja nemoteno celjenje rane.

Primeren je za suhe rane in rane z malo izločka.

KASTOR

medical dental d.o.o.

Ref.No.	Dressing Size (cm)	Pad Size (cm)	Dressings/ Box	Box/ Shipper
73230-00	5x7	2x4	10	28
73230-01	7x10	4x7	10	20
73230-02	8x15	4.5x11.5	10	20
73230-03	10x24	6x20	5	20
73230-05	10x35	6x31	5	10